Mölnlycke®

Mepilex®

WITH SAFETAC® TECHNOLOGY

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Information about materials of animal/human origin, pharmaceutical substances
Not applicable.

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
Key Accessories, Medical Devices or Non-Medical Devices, to Be Used in Combination With Medical Device Submitted For Registration.

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
Mepilex can be used under compression bandaging and in combination with gels.

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

Mepilex:

Mepilex Heel:

Mepilex can be used by lay persons under supervision of health care professionals.

1. Cleanse the wound in accordance with clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select an appropriate dressing size. Mepilex should cover the dry surrounding skin by at least 1-2 cm for the smaller sizes (sizes up to 12.5 x 12.5 cm) and 3-5 cm for the larger sizes in order to protect the surrounding skin from maceration and fixate the dressing securely. If required, Mepilex can be cut to suit various wound shapes and locations.
4. Remove the first release film and apply the adherent side to the wound.
5. Remove the remaining release film and smooth down the dressing on the skin.
Do not stretch.
6. When necessary, fixate Mepilex with a bandage or other fixation.

The dressing change interval may be several days. Change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice.

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Other information


The polyurethane foam used in the product may change colour to more yellow when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date.

If any serious incident has occurred in relation to the use of Mepilex it should be reported to Mölnlycke Health Care.


Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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

Main technical and functional characteristics

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
C – Distance grip edge release foil

D – Length of release foil overlap

E – Distance from overlapping release foil to edge

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Nom. – Nominal

Act. – Actual

Characteristics according to GOST R53498-2019	Acceptance criteria
Peeling resistance of the adhesive	

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layer, N/m: - Not less - Not more	10 1000
Functional pad sorptive capacity, cm ³ /cm ² : - not less	0.05
Functional pad wetting time, sec: - not more	10
Steam (vapor) permeability, mg/cm ² /h: - not less	1.5
Air permeability R, cm ³ /cm ² /sec: - not less	5
Atraumaticity of the functional pad (degree of adhesion to the wound surface model. % to the medical gauze): - not more	30

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Operating conditions
 Environmental conditions:
 Temperature range from +10 to +35°C.
 Relative humidity range up to 80% non-condensing (at +25°C).

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Storage and transportation conditions
 Storage conditions
 Mepilex should be stored in temperature range from 0 to +35°C.
 Mepilex should be stored in humidity range from 0 to 75%.

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Transportation conditions

Parameter	Short term/peak conditions (up to 5 hours)	Medium term conditions (up to 60 days)
Temperature [°C]	-35 to +63	0 to +36
Humidity [%]	0 to 90	0 to 75

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Package contents
 Mepilex:

Catalog number	Size, cm	Quantity in retailer	Quantity in shipping

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Mepilex Heel:

Catalog number	Size, cm	Quantity in retailer package	Quantity in shipping package
288100	13 x20	5	25

Instructions for use are published on the website <https://www.molnlycke.ru>. The link to the site is indicated on the Russian sticker.

The Russian sticker is stuck on primary (individual) and secondary (retailer) packages.

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Warranty

The manufacturer guarantees that the medical device retains its properties for the duration of the shelf life.

Disposal

The device is a potential biohazard after use and must therefore be disposed of as clinical waste.

Products coming into contact with blood and/or biological fluids are class B medical devices under medical wastes classification and may be recycled and/or disposed in accordance with applicable national regulations.

Unused products (not having contact with blood and/or biological fluids), including those with expired shelf life, belong to class A medical devices under medical waste classification and may be recycled and/or disposed in accordance with applicable national regulations.

Authorized representative of the manufacturer

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
LLC "HEALTHCARE SYSTEMS KNOWLEDGE CENTER",
127473, Moscow, Seleznevskaya st. 11A, build. 2,
floor/room/office 2/30,32/232. Tel.: +7 (499) 391-24-14;
+7 (926) 287-96-65.

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Explanation of symbols used in labelling the medical product:

Symbol	Description
	Sterilized using ethylene oxide
	Do not reuse
	CE marking
	Manufacturer
	Open here
	Caution! See instructions for use
	Keep away from moisture
	Mark waste recycling «mobius strip»
	For moderately exuding wounds

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Note! 4 first cyphers of LOT mean manufacturing date: 2 first – year (two last cyphers of year), other 2 – week of manufacturing. Ex. LOT 0952-4630 means that device was manufactured in 2009 (09=2009) on 52 week, i.e. on december.

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